

Medical Products

Quality Without Compromise











Global Gas Management Technology

Enterprises | Innovator | Medica | Westwinds

www.westernenterprises.com



Global Gas Management Technology

Our Mission

To fulfill the global demand for gas management solutions by anticipating market needs, engineering quality products & fully supporting our customers, employees and community

Our Vision

For 60 years, Western has been supplying innovative solutions in gas management systems to the

of information on Western's line of products used in the include information on available gas services, sizes, operating Additional information is available through Western's

agencies (as listed below).

to assure that safe products are delivered to the marketplace.

U.S. Food & Drug

Administration
10903 New Hampshire Avenue
Silver Springs, MD 20993
1-888-INFO-FDA

Protection Association
1 Batterymarch Park
Quincy, MA 02169
1-800-344-3555

www.westernenterprises.com

• INDUSTRIAL EQUIPMENT **PORTFOLIO**

• MANIFOLD AND SPECIALTY GAS PORTFOLIO

Gas manifold systems for industrial, medical and specialty gas applications,

• INFLATION EQUIPMENT PORTFOLIO

ISO 9001/ISO 13485 REGISTERED QUALITY SYSTEM

Disclaimer: This catalog is for informational purposes. Specific information regarding the use of these products is contained in the Operating In-



E-MAIL/INTERNET

Send orders via e-mail to customerorders@westernenterprises.com Your order will be processed.

EDI

Electronic Data Interchange (EDI) can reduce your cost of doing business. Our in-house staff will be happy to work with you to create a custom solution. Call 800-783-7890 for more information.



PHONE 800-783-7890

Press #1 to speak with a Western representative.

Please have your:

- Account number
- Part numbers, price and purchase order number
- Name of contact to answer any questions
- Special shipping instructions

BUSINESS HOURS

8:00 AM – 5:00 PM EST Monday through Friday

ORDERING

Minimum order: \$150

Freight Prepaid: \$1,250 (excluding manifolds)

	OXYGEN REGULATORS	A	OTHER RELATED EQUIPMENT OXYGEN ONLY
	OxyTOTE NG Portable Oxygen System 6-7 Compact Click-Style Oxygen Regulators 8 OctiFlo₂™ Portable Emergency Oxygen Manifold and Accessories 9 Single & Dual Stage Flow Gauge Oxygen Regulators	The state of the s	DISS Nuts, Nipples, Body Adaptors & Other
	Adjustable Pressure Regulators		OTHER GASES
	Flowmeter Regulators		DISS Nuts, Nipples & Body Adaptors 37-38 Tapered Pipe Thread Fittings
	Oxygen Cylinders and Accessories 18-19		
	Cylinder Bags and Related Equipment 20-21		GAUGES & ACCESSORIES
नीतारी नागरी से से से स्व	Cylinder Carts, Racks and Stands 22-23 FLOWMETERS		Gauges
GEN	Oxygen Flowmeters		Medical nose & Accessories49-30
л			INDEX
A CONTRACTOR OF THE PARTY OF TH	QUICK CONNECTSOxygen Quick Connects27Medical Air Quick Connects28Nitrous Oxide Quick Connects29Suction Quick Connects30		CGA Gas Assignments Chart



OxyTOTE NG System Advantages

Western designed the Next Generation Integrated Valve Regulator, the OxyTOTE NG, to meet the demands of both EMS emergency services as well as in-hospital patient transport. It is small, lightweight, and ergonomically designed. Every aspect of its intended use was carefully considered during the design process, including its length, which allows the OxyTOTE NG to fit into most EMS duffle bags as well as horizontally on most hospital gurneys.

The OxyTOTE NG's SureClick $^{\text{TM}}$ feature locks in oxygen flow with accuracy. The high impact shroud protects the regulator, gauge, and hose barb from accidental damage.

Features

- Regulator and cylinder valve integrated into one unit
- User can initiate flow in two easy steps
- Compact design with all brass body
- Built-in Safety Burst Disk prevents over pressurization
- 1¹/₂" diameter oxygen content gauge continuously indicates cylinder pressure
- CGA-540 Fill Port with Safety Check Valve
- Filtration throughout for maximum safety and reduced particle contamination
- Protective shroud ensures safety and lasting performance
- 50 PSI DISS Demand Valve Outlet
- 5 yr limited warranty

OxyTOTE NG Regulator Specifications

(OxyTOTE NG Regulator only)

Size	5.240" x 4.375" (5.989" x 4.375" including threads)
Weight	2.14 lbs.
OxyTOTE NG Flow*(LPM)	0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15, 25

^{*}The 80X Series excludes 8 LPM flow setting.



OxyTOTE NG System Advantages

Part#	Description
MNR-600	OxyTOTE NG regulator head only, with DISS Demand Outlet, 3/4"–16UNF thread (aluminum cylinder)1
MNR-800	OxyTOTE NG regulator head only, with DISS Demand Outlet, 3/4"–16UNF thread (aluminum cylinder). No shut off valve ¹
MNS-602	OxyTOTE NG system includes D size aluminum cylinder ¹
MNS-603*	OxyTOTE NG system includes E size aluminum cylinder ¹
MNS-802	OxyTOTE NG system includes D size aluminum cylinder. No shut off valve ¹
MNS-803*	OxyTOTE NG system includes E size aluminum cylinder. No shut off valve ¹
PX-1044	Humidification Attachment Sleeve
PX-1045	Humidification Hose

^{*}MR Conditional for use 1.5 ft. or more outside of the bore of an MR machine under the following conditions: Static magnetic field of 3.0 Tesla or less and a spatial gradient field of 720 gauss/cm or less.

OxyTOTE NG PX-1044
with optional
humidification
attachment.

OxyTOTE NG - The Next Generation of Portable Oxygen Delivery Systems

Finally, an oxygen delivery system that operates as easy and quickly as you do. OxyTOTE NG delivers oxygen in two easy steps. No hassles with wrenches, washers or attachments – just a simple way for you to do your job and ensure your patient is getting the care they need quickly and efficiently.

Confidence – Designed and manufactured within a quality system registered to both ISO9001 and ISO13485 ensures OxyTOTE NG's performance, longevity, and total satisfaction.

Made in USA – OxyTOTE NG is manufactured and assembled in Westlake, OH.







OXYTOTE SYSTEM SPECIFICATIONS

SPECIFICATIONS	MN5-002	MN2-003
Cylinder Size	425 liter "D" Aluminum	679 liter "E" Aluminum
System Size	21-5/8"H x 4-1/2"W	30-5/8"H x 4-1/2"W
System Weight	7.4 lbs.	9.8 lbs.
Regulator Flow (SLPM)	1/2 - 25	1/2 - 25

¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Designed for pediatric, homecare, hospital and emergency applications where durability and accuracy matter.

OPA SERIES FEATURES:

- Accurate Click-Style flow control
- Maximum inlet pressure 3000 PSI
- Sintered filter for additional safety and to extend regulator life
- Internal reseating relief valve protects against over-pressurization
- 1-1/2" diameter gauge
- Tapered hose barb outlet
- 2 year limited warranty
- Durable compact aluminum body



Compact Click-Style Regulator with CGA-870 Yoke Inlet



OPA-520 Compact Click-Style Regulator with handtight CGA-540 Nut and Nipple Inlet

OPA SERIES COMPACT REGULATORS

Part #	Flow Settings		
CGA-870 Yok	CGA-870 Yoke		
OPA-810	1/64, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1-1/2, 2 LPM ¹		
OPA-820	1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15 LPM ¹		
OPA-840	1/2, 1, 1-1/2, 2, 2-1/2, 3, 4, 5, 6, 8 LPM ¹		
OPA-850	1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 LPM ¹		
CGA-540 Nut and Nipple – Handtight			
OPA-510	1/64, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1-1/2, 2 LPM ¹		
OPA-520	1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15 LPM ¹		
OPA-540	1/2, 1, 1-1/2, 2, 2-1/2, 3, 4, 5, 6, 8 LPM ¹		

EXTENDED CLICK-STYLE REGULATORS

OPL SERIES FEATURES:

- Same as OPA Series PLUS:
- Extended body with DISS 1240 outlet for use with a humidifier bottle

OPL SERIES EXTENDED REGULATORS

Part #	rt # Flow Settings Outlet Cor	
CGA-870 Yoke		
OPL-821	1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15 LPM ¹	DISS 1240



Extended Click-Style Regulator with DISS 1240 Outlet Connection (shown with humidifier bottle) not available from Western

WESTERN MINI OXYGEN REGULATORS

Western Mini Oxygen Regulators are everything you need for homecare in the smallest, lightest, most affordable regulator we've ever made!

FEATURES:

- Click-Style flow control for ease of use
- Compact and lightweight aluminum body to aid portability
- Brass high-pressure insert to assure safety
- Sintered filter for additional safety and to extend regulator life
- 1-1/2" diameter gauge is easy to read
- Maximum inlet pressure 3000 PSI
- Available in CGA-870 yoke style only
- 1 year limited warranty

Part #	Flow Settings	Outlet Connection
MLR-115	1/4, 1/2, 1, 1-1/2, 2, 3, 4, 5, 6, 7, 8, 10, 15 ¹	1/4" OD hose barb



¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

OCTI-FLO₂ Portable Emergency Oxygen Manifold



OCTI-FLO₂ Portable Emergency Oxygen Manifold

- Used for mass casualty, EMS, and triage situations
- Oxygen delivery for up to 8 patients
- Flowmeters are individually controlled
- 8 easy-to-read Click-Style flowmeters from 0-25 lpm
- Pressure gauge on manifold
- Standard DISS 1240 inlet and a 100 PSI relief valve safety feature
- Lightweight and compact durable aluminum construction
- 4 Cannulas, 8 breathing masks, and compact carrying case included
- Affordable high quality craftsmanship for years of dependable service
- 2 year limited warranty

Optional Components to Complete the Emergency Oxygen Delivery System |

- 0-100 PSI Regulator
- DISS 20' Oxygen Hose



Accessories for Other EMS Applications

- Wrenches • T-handles
- Replacement Parts
- Post Valves
- Washers

MCW - 2BC



¹▲**WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

SINGLE & TWO-STAGE OXYGEN FLOW GAUGE REGULATORS

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauges
- 6 year limited warranty

COMPACT MODELS (MR SERIES)

- Single-stage with CGA-870 swivel yoke to add convenience for orientation
- 1-1/2" diameter gauges

TWO STAGE MODELS:

• Two stage design for applications where flow control accuracy is critical such as calibration and testing



HELPFUL TIPS

First stage reduces high pressure to 170 PSI; second stage reduces to provide greater accuracy at lower flow settings

MODEL #	FLOW RANGE	INLET	
SINGLE-STAGE			
M1-870-5FG	1/2-5 LPM	Yoke ¹	
M1-540-8FG	1-8 LPM	Nut and Nipple 1	
M1-540-15FG	2-15 LPM	Nut and Nipple ¹	
M1-870-15FG	2-15 LPM	Yoke ¹	
MR-870-15FG	2-15 LPM	Yoke ¹	
PEDIATRIC SINGLE STAGE			
M1-540-PED3	1/32- 1/2 LPM	Nut and Nipple 1	
M1-870-PED3	1/32- 1/2 LPM	Yoke ¹	
M1-540-PED	1/8-1 LPM	Nut and Nipple 1	
M1-870-PED	1/8-1 LPM	Yoke ¹	
MR-870-PED	1/8-1 LPM	Swivel Yoke ¹	
TWO-STAGE			
M2-540-15FG	2-15 LPM	Nut and Nipple ¹	



M1-870-15FG

Single-Stage Regulator with CGA-870 Connection



Two-Stage Regulator with CGA-540 Connection



¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



M1-910-PG Single-Stage Adjustable Regulator with CGA-910 Yoke



M1-326-PGH Single-Stage Adjustable Regulator with CGA-326 Handtight Nut



M1-500-PG

FEATURES:

- Chrome-plated brass with all brass high-pressure chamber
- Maximum inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauges
- Dual scale: PSI and kPa
- Outlets are DISS gas specific
- 6 year limited warranty

ADJUSTABLE SINGLE-STAGE MODELS

GAS/MODEL	DELIVERY PRESSURE	CGA INLET CONNECTION
CARBON DIOXIDE		
M1-320-PG	0 - 100 PSI	320-Nut and Nipple ¹
M1-320-PGH	0 - 100 PSI	320-Handtight Nut and Nipple ¹
M1-940-PG	0 - 100 PSI	940-Yoke ¹
NITROUS OXIDE		
M1-326-PG*	0 -100 PSI	326-Nut and Nipple ¹
M1-326-PGH*	0 -100 PSI	326-Handtight Nut and Nipple 1
M1-910-PG*	0 -100 PSI	910-Yoke ¹
MEDICAL AIR		
M1-346-PG	0 -100 PSI	346-Nut and Nipple ¹
M1-346-PGH	0 -100 PSI	346-Handtight Nut and Nipple 1
M1-950-PG	0 -100 PSI	950-Yoke ¹
OXYGEN		
M1-540-PG	0 -100 PSI	540-Nut and Nipple ¹
M1-540-PGH	0 -100 PSI	540-Handtight Nut and Nipple ¹
M1-870-PG	0 -100 PSI	870-Yoke ¹
NITROGEN		
M1-580-PG	0 -100 PSI	580-Nut and Nipple 1
M1-960-PG	0 -100 PSI	960-Yoke ¹
HELIUM		
M1-930-PG	0 -100 PSI	930-Yoke ¹
MIXTURES		
M1-500-PG	0 -100 PSI	500-Nut and Nipple ¹
M1-500-PGH	0 -100 PSI	500-Handtight Nut and Nipple ¹

^{*1000} PSI gauge color – coded blue See full CGA chart in back of catalog.

¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauges
- 50 PSI preset output
- 6 year limited warranty

M1A

PRESET WITH WITH 90° HORIZONTAL OUTLET OUTLET

GAS/MODEL OXYGEN		CGA	CGA INLET
OXYGEN		COAL	CGA INLEI
OXIOLIV			
M1A-540-P	M1-540-P	540	Nut and Nipple 1
M1A-540-PH	M1-540-PH	540	Handtight Nut and Nipple 1
M1A-870-P	M1-870-P	870	Yoke ¹
N/A	MR-870-P	870	Yoke (Compact Model) 1
CARBON DIOXID	E		
	M1-320-P	320	Nut and Nipple 1
M1A-940-P	M1-940-P	940	Yoke ¹
NITROUS OXIDE			
	M1-326-P	326	Nut and Nipple 1
M1A-910-P	M1-910-P	910	Yoke ¹
MEDICAL AIR			
M1A-346-P	M1-346-P	346	Nut and Nipple 1
M1A-346-PH	M1-346-PH	346	Handtight Nut and Nipple 1
M1A-950-P	M1-950-P	950	Yoke ¹
NITROGEN			
	M1-580-PH	580	Handtight Nut and Nipple 1
	M1-580-P	580	Nut and Nipple 1
	M1-960-P	960	Yoke ¹
MEDICAL GAS M	IXTURES		
	M1-280-P	280	Nut and Nipple ¹
	M1-500-P	500	Nut and Nipple 1
	M1-500-PH	500	Handtight Nut and Nipple 1
	M1-890-P	890 Yoke ¹	
MEDICAL HELIUM			
	M1-580A-PH	580	Handtight Nut and Nipple 1

See full CGA chart on inside back cover.



M1A-540-P

Single-Stage Preset Regulator with CGA-540 Nut & Nipple Inlet with Horizontal Outlet



M1-540-P

Single-Stage Preset Regulator with CGA-540 Nut & Nipple Inlet with 90° Outlet



M1A-870-P

Single-Stage Preset Regulator with CGA-870 Yoke Connection with Horizontal Outlet

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



Single-Stage Preset Flowmeter Regulator with CGA-540 Nut and Nipple Inlet



Single-Stage Flowmeter Regulator with CGA-870 Yoke

OXYGEN MR CONDITIONAL SLIM-LINE FLOWMETER

MODEL	FLOW RANGE	REGULATOR PRESSURE	INLET
NMR-870-15FM	1/2 - 15 LPM	50 PSI preset	CGA-870 Yoke ¹

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum rated inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauge
- 15 LPM Oxygen models provide 75 LPM flood minimum
- 7 LPM models restricted to 10 LPM maximum
- Polycarbonate inner and outer tubes provide greater accuracy and durability
- Includes plastic swivel hose barb (not shown)
- 6 year limited warranty

SINGLE STAGE MODELS

GAS/MODEL	FLOW RANGE	INLET CONNECTION
OXYGEN		
M1-540-8FM	1/2-7 LPM	CGA-540 Nut and Nipple ¹
M1-540-8FMH	1/2-7 LPM	CGA-540 Handtight Nut and Nipple 1
M1-870-8FM	1/2-7 LPM	CGA-870 Yoke ¹
M1-540-15FM	1/2 - 15 LPM	CGA-540 Nut and Nipple ¹
M1-540-15FMH	1/2-15 LPM	CGA-540 Handtight Nut and Nipple 1
M1-870-15FM	1/2 - 15 LPM	CGA-870 Yoke ¹
M1-540-15FMPTO	1-15 LPM	CGA-540 Nut and Nipple w/Power Take Off ¹
MEDICAL AIR		
M1-346-15FM	1-15 LPM	CGA-346 Nut and Nipple ¹
M1-346-15FMH	1-15 LPM	CGA-346 Handtight Nut and Nipple 1
M1-950-15FM	1-15 LPM	CGA-950 Yoke ¹
CARBON DIOXIDE		
M1-320-12FM	1/2 - 12 LPM	CGA-320 Nut and Nipple ¹
M1-320-12FMH	1/2 - 12 LPM	CGA-320 Handtight Nut and Nipple 1
M1-940-12FM	1/2 - 12 LPM	CGA-940 Yoke ¹

See full CGA chart on inside back cover

FLOWMETER FEATURES:

- Polycarbonate inner and outer tubes
- Calibrated at 50 PSI inlet
- DISS 1240 outlet

REGULATOR FEATURES:

- All-Brass body and high-pressure chamber
- Internal reseating relief valve protects against over-pressurization
- 2" diameter gauge
- Sintered filter for additional safety and extended regulator life
- Durable neoprene diaphragm
- 3000 PSI maximum inlet pressure
- DISS 1240 outlet
- 50 PSI preset outlet

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauges
- OUTLET: Needle Valve with 1/8" Hose Barb (all models)
- 6 year limited warranty

PART #	CGA	INLET
M1-280-PGB	280	Nut and Nipple 1
M1-320-PGB	320	Nut and Nipple 1
M1-950-PGB*	950	Yoke ¹
M1-500-PGB	500	Nut and Nipple 1
M1-973-PGB*	973	Yoke ¹
M1-540-PGB	540	Nut and Nipple 1
M1-870-PGB*	870	Yoke ¹

Note: For Handtight Nut add "H" to item number.

See full CGA chart on inside back cover





Regulator with CGA-870 Yoke Inlet

^{*}Handtight option not available

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum inlet pressure 3000 PSI
- Durable Neoprene diaphragm
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauges
- Brass outlet shut-off valve
- External relief valve
- 1/4" male NPT outlet

HIGH PURITY MODELS

MODEL	DELIVERY RANGE	MAX RATED FLOW SCFH AIR AT 2000 PSI*	INLET CONNECTION
MSH15XXX	0 - 15 PSI	500 SCFH	Nut and Nipple or Yoke ¹
MSH45XXX	0 - 45 PSI	650 SCFH	Nut and Nipple or Yoke ¹
MSH80XXX	0 - 80 PSI	950 SCFH	Nut and Nipple or Yoke ¹
MSH180XXX	0 - 180 PSI	1000 SCFH	Nut and Nipple or Yoke ¹
MSH450XXX	0 - 450 PSI	1000 SCFH	Nut and Nipple or Yoke ¹
MSG80XXX	0 - 80 PSI	2050 SCFH	Nut and Nipple or Yoke ¹
MSG180XXX	0 - 180 PSI	2750 SCFH	Nut and Nipple or Yoke ¹

^{*}Maximum rated flow per CGA E-4 1994 at the stated inlet pressure



MSH180346

Single-Stage High Purity Medical Air Regulator with CGA-346 Nut and Nipple Inlet

ORDERING EXAMPLE

Order by replacing X's in model number with desired CGA number from the Gas/Inlet Connection chart. Example: MSH45XXX converts to MSH45320 for a Carbon Dioxide regulator with a CGA-320 nut and nipple inlet connection.

GAS/INLET CONNECTIONS

NUT & NIPPLE CGA	YOKE CGA	GAS ASSIGNMENT	
280	880	Carbon Dioxide and Oxygen mixtures (CO ₂ not over 7%) ¹	
280	890	Helium and Oxygen mixtures (He not over 80%) ¹	
280	890	Nitrogen and Oxygen mixtures (O ₂ over 23.5%) ¹	
280	890	Xenon and Oxygen (O ₂ over 20%) ¹	
320	940	Carbon Dioxide ¹	
326	910	Nitrous Oxide ¹	
346	950	Medical Air ¹	
500	940	Carbon Dioxide and Oxygen mixtures (CO ₂ over 7%) ¹	
500	930	Helium and Oxygen mixtures (He not over 80%) ¹	
500	973	Carbon Dioxide, Oxygen and Nitrogen ¹	
500	973	Lung Diffusion Mixtures ¹	
500	973	Nonflammable, Noncorrosive Diagnostic and Medically Related Gas Mixtures ¹	
540	870	Oxygen ¹	
580	960	Nitrogen ¹	
580	930	Helium ¹	

Not all gases are available in various delivery pressure ranges.

EASY GUIDE TO WESTERN MEDICA PART NUMBERS

MEDICAL GAS FITTINGS

HOW IT WORKS:

M24-1 (M=Medical 24=DISS 1240 1= Part#)

(M=Medical 54=CGA-540 1=Part #) Medical Oxygen Cylinder Nut

(M=Medical 87=CGA-870 1=Part #)

Medical Oxygen Yoke

MEDICAL GAS REGULATORS

Part Number-M1-540-FG

Medical Oxygen DISS Nut

(M=Medical 1=Single Stage 540=CGA-540 FG=Flow Gauge)
Description: Medical Single-Stage Oxygen Flow Gauge Regulator

¹▲**WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

BREATHING MIXTURE AND MEDICAL MIXTURE REGULATORS

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauges
- 6 year limited warranty





BREATHING MIXTURES AND MEDICAL MIXTURES

		SINGLE STAGE ADJUSTABLE REGULATORS 0-100 PSI				
CGA		PART #	INLET	DISS OUTLET	OUTLET PART #	
	BREATHING MIXTURES					
280	CO ₂ & O ₂ (CO ₂ not over 7%)	M1-280-PG	Nut and Nipple ¹	1020	M02-5	
880		M1-880-PG	CGA-880 Yoke ¹	1200A	M20-5	
890	He & O ₂ (He not over 80%)	M1-890-PG	CGA-890 Yoke ¹	1240	M24-5	
890	N ₂ & O ₂ (O ₂ over 23.5%)	M1-890-PG	CGA-890 Yoke ¹	1240	M24-5	
890	Xe & O ₂ (O ₂ over 20%)	M1-890-PG	CGA-890 Yoke ¹	1240	M24-5	
	MEDICAL MIXTURES					
940	CO ₂ & O ₂ mixtures (CO ₂ , over 7%)	M1-940-PG	CGA-940 Yoke ¹	1080A	M08-5	
930	He & O ₂ (He over 80%)	M1-930-PG	CGA-930 Yoke ¹	1060A	M06-5	
500	CO ₂ & O ₂ & N ₂	M1-500-PG	Nut and Nipple ¹	1020	M02-5	
973	LUNG DIFFUSION MIXTURES	M1-973-PG	CGA-973 Yoke ¹	1020	M02-5	
	BLOOD GAS MIXTURES					
500	Nonflammable, Noncorrosive Gas	M1-500-PG	Nut and Nipple ¹	1020	M02-5	
973	Mixtures Labeled as Drugs or Medical Devices and not having another connection assignment	M1-973-PG	CGA-973 Yoke ¹	1020	M02-5	

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

MAGNETIC RESONANCE (MR) CONDITIONAL* REGULATORS



Slim Line-Flowmeter



Single-Stage Flowmeter Regulator

FLOWMETER FEATURES:

- Slim-Line Flowmeter assembled with non-magnetic components
- Anodized aluminum body with brass adjustment knob
- Polycarbonate inner and outer tubes
- Calibrated at 50 PSI inlet
- DISS 1240 outlet

REGULATOR FEATURES:

- Flowmeter Regulator assembled with non-magnetic components
- Maximum inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauge
- All-Brass body and high-pressure chamber
- DISS 1240 outlet
- 50 PSI outlet

FLOWMETER MODELS

GAS	MODEL	FLOW RANGE	INLET PRESSURE	INLET FEMALE CONNECTION
Oxygen	FMR100	1/2 - 15 LPM	50 PSI	1/8 - 27 NPT ¹
Oxygen	FMR103	1/2 - 15 LPM	50 PSI	DISS 1240 Nut And Nipple ¹

FLOWMETER REGULATOR MODELS

G	AS M	ODEL FLOW RA	ANGE INLET PRES	SURE INLET CONNECTION
Ox	vaen NMR-8	370-15FM 1/2 - 15	LPM 50 PSI	CGA - 870 Yoke ¹

^{*}The products on this page are MR conditional for use outside the bore of an MR machine under the following conditions: Static magnetic field of 3.0 Tesla or less; Spatial gradient field of 3.3 Tesla/meter or less.



¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CYLINDER VALVES (3000 PSI)

3/4"-16 UNF-2A



Note: The cylinder valves listed are configured with proper relief mechanisms for use at 2015 and 2216 PSI. Do not use where cylinder pressure ratings exceed these pressures.

ALUMINUM CYLINDER FEATURES:

- Aluminum cylinders are 40% lighter than steel
- Cylinders are cleaned for medical oxygen service
- Cylinder finish is easy to clean and maintain
- Straight thread design minimizes thread damage and with o-ring provides positive seal
- Order quantity minimum of 6

POST VALVE SEALS

Post Valve Seal Holder¹ **MPV-4027** (Plastic Clip Dust Cover) see photo above

Note: See page 34 for more oxygen accessories



M6 (B)

1	Portable Aluminum Cylinders with CGA-870 Post Valves*	11-3/4" x 3" • 2.2 lbs. 164 liters • 2216 PSI 144/Pallet • 5.8 cu ft
		WPX-8700-1W WPX-8700-1T

^{*}Weights listed are without valves. Add 1/2 lb for each valve.

¹▲**WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

2 WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

18

CYLINDER WRENCHES



Size shown not proportionally correct.

WPX-8703-1T

Spigot Valve

Spigot Valve

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

SHOULDER PACKS

Part #	Description
PX-1011	Black Nylon Bag for C Cylinder
PX-1013	Black Nylon Bag for D Cylinder
PX-1014	Black Nylon Bag for E Cylinder
PX-1010	Horizontal Carrying Bag (duffel style), vented Holds M4 (A), M6 (B), M6A (ML6), and M9 (C) Cylinders







OTHER EQUIPMENT

Part #		
MW-87	Halar Yoke Washer (one time use only)	
2544-2	Brass and Viton Yoke Washer ¹	
PX-1850	Dual Lumen Cannula, 5'	
PX-1000	Single Lumen Cannula, 7' Call to order a variety of pressure gauges. See page 46	
PX-1002	Adult O ₂ Mask w/ tubing	
M24-45P	Nut and Nipple to Tapered Barb 50/bag (green)	
M24-45PGR	Same as above 50/bag (Grey)	
M24-45PW	Same as above 50/bag (white)	
M24-45PY	Same as above 50/bag (yellow)	
M00-3	T-Handle ¹	
M00-9	T-Handle with 6" grip ¹	
MTT-1K	1/4" PTFE Tape	
MTT-2K	1/2" PTFE Tape	
MTT-3K	Formula 8 PTFE Paste	





2544-2





PX-1850

PX-1002



M24-45PW









¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CYLINDER CARTS²



2▲**WARNING:** This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CYLINDER RACKS AND STANDS 2

Part #	Cylinder	Capacity	Width (in.)	Depth (in.)	Height (in.)	Weight (lbs.)
PX-1025	C,D or E	6	16	10.5	7.25	10
PX-1027	C,D or E	12	20.5	15.8	7	13.18
PX-1034	D or E	40	40	25	20	75
MCS	C,D or E	1	9.4	9.4	4	2
HCS2-M	M60, M or H	1	19	18.25	11	14
HCS3-M	M60 or M	1	19	18.25	11	5
















MCS

HCS2-M

²▲WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

INLET	FLOWMETERS WITH FITTING
Flowmeter with Standard Inlet (1/8 NPT Female) ¹	FM101 1/2–15 LPM FM201 1/2–7 LPM
1/4 NPT Male ¹	FM102 1/2–15 LPM FM202 1/2–7 LPM
DISS Hex Nut and Nipple ¹	FM103 1/2–15 LPM FM203 1/2–7 LPM
DISS Handtight Nut and Nipple*	FM104 1/2–15 LPM FM204 1/2–7 LPM

 $^{^{\}star}$ Canadian color-coding available on most models, add $^{\prime\prime}$ C $^{\prime\prime}$ to part number.

FLOWMETER FEATURES:

- Full line of inlet connections
- Optional Power Take-Off units provide 50 PSI additional outlet pressure
- Polycarbonate inner and outer tubes provide greater durability
- Calibrated at 50 PSI inlet
- Double units for multi-patient use
- DISS 1240 outlet
- 6 year limited warranty

OPTIONS:

• Order 2052 for safety chain assembly.



FM101

¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

FLOWMETER	S WITH FITTING	INLET
FM106	1/2–15 LPM ¹	Oxequip® Adaptor
FM107	1/2–15 LPM ¹	Ohmeda® Adaptor*
	1/2–15 LPM ¹ 1/2–7 LPM ¹	Puritan® Adaptor
FM109	1/2–15 LPM ¹	Chemetron® Adaptor*

^{*}Canadian color-coding available on most models, add "C" to part number. Example FM107C

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

FM - Brass Body

INLET	FLOWMETERS WITH FIT	TING
Flowmeter Without Fitting Standard Inlet	MEDICAL AIR 1-15 LPM	CARBON DIOXIDE 1/2-12 LPM
1/8 NPT Female ¹	FM601	FM801
1/4 NPT Male ¹	FM602	FM802
DISS Hex Nut and Nipple ¹	FM603	FM803
DISS Handtight Nut and Nipple*	FM604	NA
Ohmeda® Adaptor*1	FM607	NA
Chemetron® Adaptor*1	FM609	NA

^{*}Canadian color-coding available on most models, add "C" to part number. Example FM604C

FLOWMETER FEATURES:

- Full line of inlet connections
- Optional Power Take-Off (PTO) units provide 50 PSI additional outlet pressure
- Polycarbonate inner and outer tubes provide greater durability
- Calibrated at 50 PSI inlet
- Double units for multi-patient use
- 6 year limited warranty


OPTIONS:

• Order 2052 for safety chain assembly.

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Chemetron® Male Adaptor*	MQ101 MQ101-1 MQ101-2 MQ101-3 MQ101-4 MQ101-5 MQ101-6	1/8 NPT Male ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Female ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Chemetron® Female Adaptor*	OR101 OR101-1 OR101-2 OR101-3 OR101-4 OR101-5 OR101-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Male Adaptor*	MQ201 MQ201-1 MQ201-2 MQ201-3 MQ201-4 MQ201-5 MQ201-6	1/8 NPT Male ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Female ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Female Adaptor*	OR201 OR201-1 OR201-2 OR201-3 OR201-4 OR201-5 OR201-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Oxequip® Male Adaptor	MQ301	1/8 NPT Male ¹
Puritan® Male Adaptor	MQ401	1/8 NPT Male ¹

50 PSI Nominal Working Pressure










^{*}Canadian color–coding available on most models. Add "C" to part number. Example: MQ101-4C

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



DISS Nut and Nipple

Chemetron® Male Style*	MQ102 MQ102-1 MQ102-2 MQ102-3 MQ102-4 MQ102-5 MQ102-6	1/8 NPT Male ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Female ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Chemetron® Female Style*	OR102 OR102-1 OR102-2 OR102-3 OR102-4 OR102-5 OR102-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Male Style*	MQ202 MQ202-1 MQ202-2 MQ202-3 MQ202-4 MQ202-5 MQ202-6	1/8 NPT Male ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Female ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Female Style*	OR202 OR202-1 OR202-2 OR202-3 OR202-4 OR202-5 OR202-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Oxequip® Male Style	MQ302	1/8 NPT Male ¹
Puritan® Male Style	MQ402	1/8 NPT Male ¹

50 PSI Nominal Working Pressure

^{*}Canadian color-coding available on most models. Add "C" to part number. Example: MQ204-2C

¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Chemetron® Female Style*	OR103 OR103-1 OR103-2 OR103-3 OR103-4 OR103-5 OR103-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Male Style*	MQ203 MQ203-1 MQ203-2 MQ203-3 MQ203-4 MQ203-5 MQ203-6	1/8 NPT Male ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Female ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Female Style*	OR203 OR203-1 OR203-2 OR203-3 OR203-4 OR203-5 OR203-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Oxequip® Male Style	MQ303	1/8 NPT Male ¹
Puritan® Male Style	MQ403	1/8 NPT Male ¹



^{*}Canadian color–coding available on most models. Add "C" to part number. Example: OR203-2C



Male



1/4" Hose Barb Female



DISS Male









¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



3	
	1/4" NPT Male





DISS Nut and Nipple

Chemetron® Male Style*	MQ104 MQ104-1 MQ104-3	1/8 NPT Male ¹ 1/4 Hose Barb ¹ 1/4 NPT Male ¹
Chemetron® Female Style*	OR104-1 OR104-2 OR104-3 OR104-4 OR104-5 OR104-6	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Male ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Male Style*	MQ204 MQ204-1 MQ204-2 MQ204-3 MQ204-4 MQ204-5 MQ204-6	1/8 NPT Male ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹ 1/8 NPT Female ¹ Handtight Nut and Nipple ¹ DISS Nut and Nipple ¹
Ohmeda® Female Style*	OR204 OR204-1 OR204-2 OR204-3	1/8 NPT Female ¹ 1/4 Hose Barb ¹ DISS Male ¹ 1/4 NPT Male ¹
Oxequip® Male Style	MQ304	1/8 NPT Male ¹
Puritan® Male Style	MQ404	1/8 NPT Male ¹

50 PSI Nominal Working Pressure

*Canadian color–coding available on most models. Add "C" to part number.

¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

DISS 1240	200 PSI Max.				
NUTS (9/1	6-18 Thread)	NIPPLES (DISS 1240)	BODY AD	APTORS (DISS 1240)
M24-1	11/16" Hex ¹	M24-10	1/8" NPT Male x 1-1/2" w/Wrench Flats ¹	M24-4	1/8" NPT Male ¹
M24-2	Knurled Handtight ¹	M24-11	1/8" NPT Male x 1-1/2" w/Wrench Flats and O-Ring ¹	M24-5	1/4" NPT Male ¹
M24-2P	Handtight	M24-12	Tapered 3/16" to 3/8" Hose Barb ¹	M24-6	1/8" NPT Female ¹
	w/Green Collar ¹	M24-13	Tapered 3/16" to 3/8" Hose Barb w/O-Ring ¹	M24-7	1/4" NPT Female ¹
M24-2PA	Handtight Nut Ass'y 1/8" NPT Green ¹	M24-14	1/4" Hose Barb ¹	M24-8	1/4" Hose Barb ¹
M24-2PH	Hex Handtight	M24-15	1/4" Hose Barb w/O-Ring ¹	M24-9	3/16" Hose Barb ¹
MZ-7 ZI II	w/Green Collar ¹	M24-16	3/16" Hose Barb ¹		
M24-2PW	Handtight Canada	M24-18	O-Ring Replacement ¹		
	w/White Collar ¹	M24-41	1/8" NPT Male x 1" Broached for Allen Wrench ¹		
M24-33	DISS 1240F x 1/4" NPT Male ¹	M24-42	1/8" NPT Male x 1" Broached w/O-Ring ¹		
		M00-6	Retaining Ring ² M24-10	200	M24-5
M24-2	M24-1		M24-13	auc	M24-7
M24-2P	M24-2PA		M24-14		M24-8

Important: Hose used on straight barbed hose fittings should be secured with crimp-style hose ferrules or proper clamps. Nipples with O-Ring should be used when using Handtight nuts to insure proper seal.

OXYGEN DISS SPECIAL ADAPTORS

DISS 1240 9/16" Thread 200 PSI	Max.			
90°	STRAIGHT		MISCELLANEOU	IS
M24-27 Female to Male ¹	M24-32 Fe	emale to 1/8" NPT Male ¹	M24-24 Fema	e to 1/4" Hose Barb ¹
M24-28 Female to 1/8" NPT Male ¹	M24-33 Fe	emale to 1/4" NPT Male ¹	M24-25 Fema	e to Female ¹
M24-29 Female to 1/4" NPT Male ¹	M24-34 Fe	emale to 1/8" NPT Female ¹	M24-26 Male	to Male ¹
M24-30 Female to 1/8" NPT Female ¹	M24-35 Fe	emale to 1/4" NPT Female ¹	M24-47 Femal	e Handtight to Tapered ¹
M24-31 Female to 1/4" NPT Female ¹				
M24-39 Male to 1/4" NPT Male ¹				
M24-43 Male to Female, with 1-3/4" Nipple Extension	1			
M24-27	<u> </u>	M24-36		M24-25
M24-29		M24-33		M24-26
M24-39	(M24-35		M24-24

¹ **_____WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

OXYGEN DUST CAPS, GAS TIGHT CAPS, SWIVEL CONNECTIONS

DISS 1240	9/16-18 UNF Thread	200 PSI Max.		
DUST CAPS	(NOT FOR GAS SEAL)	GAS TIGHT CAPS	PLASTIC SWI	VEL CONNECTORS*
M24-48H	Black Cap w/Chain and Ring ²	M24-20 Chrome Cap w/Chain and Ring ²	M24-45P	Nut and Nipple to Tapered Barb (Green)*
M24-485	Black Cap w/Chain and Eyelet ²		M24-45PGR	Nut and Nipple to Tapered Barb (Grey)*
			M24-45PW	Nut and Nipple to Tapered Barb (White)*
			M24-45PY	Nut and Nipple to Tapered Barb (Yellow)*
	M24-48H	M24-20		M24-45P M24-45PW

^{*}Sold in bags of 50



ANODIZED ALUMINUM MULTIPLE PORT BLOCK (oxygen only) Two 1/4 NPT Outlets 200 PSI Max.



M00-70

Anodized Aluminum 200 PSI

M00-70 1/8" NPT Female Inlet to 1/4" NPT Female Outlets²
M00-71 1/4" NPT Female Inlet to 1/4" NPT Female Outlets²

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

OXYGEN	DISS 1240 9/16" Thr	ead 200	PSI Max		
DEMAND VA	ALVES	CHECK VA	LVES	REVERSE	FLOW VALVES
500 SCFH Flo	w @ 50 PSI	50 PSI		1000 SCFH	H flow @ 50 PSI
M24-40	1/8" NPT Male to DISS ¹	M24-4CV	1/8" NPT Male 50 PSI ¹	MCV-1	Flow into 9/16-18 Male Inlet to Female Outlet ¹
M24-50	1/4" NPT Male to DISS ¹	M24-5CV	1/4" NPT Male 50 PSI ¹	MCV-2	Flow into 9/16-18 Male inlet to Male Outlet ¹
		MCV-3	Flow into 1/4 NPT Male Inlet to 9/16-18 Male Outlet ¹		
Cimin 5	M24-50	•	M24-5CV		MCV-1 MCV-2
					MCV-3
Adaptors					
MD-1	1/8" NPT Male to 1/8" NI	PT Female ¹			Must be used with demand valve when plunger
MD-2	1/4" NPT Male to 1/8" NI	PT Female ¹		hole depth	is not sufficient to clear forward travel of plunger ection is made. Always check for full flow after
MD-3	1/4" NPT Male to 1/4" NI	PT Female ¹	MD-3	installation.	

Demand Valves

Gas is checked w/out connection being made. Flow occurs when proper connection is made.

Check Valves

Gas flows into and through the DISS portion of the fitting and is checked in the opposite direction. Does not operate as a demand valve.

Reverse Check Valves

For inline gas flow. Prevents backflow of gases.

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



BOTTOM P	ORT				BACK PORT	Г			
MG-C1	4000 PSI	1-1/2" OD	1/8" NPT	Chrome ¹	MG-C6*	15 LPM	1-1/2" OD	1/8" NPT	Chrome ¹
MG-C2	4000 PSI	2" OD	1/4" NPT	Chrome ¹	MG-C23	4000 PSI	1-1/2" OD	1/8" NPT	Chrome ¹
MG-C3	100 PSI	1-1/2" OD	1/8" NPT	Chrome ¹	MG-C67*	1-8 LPM	1-1/2" OD	1/8" NPT	Chrome ¹
MG-C4	100 PSI	2" OD	1/4" NPT	Chrome ¹	MG-C75	3000 PSI	1-1/2" OD	1/8" NPT	Chrome ¹
MG-C5	15 LPM*	1-1/2" OD	1/8" NPT	Chrome ¹	MJ-C1P	3000 PSI	1-1/2" OD	1/8" NPT	Stainless ²
MG-C7	15 LPM*	2" OD	1/4" NPT	Chrome ¹	MJ-C2	4000 PSI	1-1/2" OD	1/8" NPT	Stainless
MG-C8	5 LPM*	2" OD	1/4" NPT	Chrome ¹					
MG-C31	8 LPM*	2" OD	1/4" NPT	Chrome ¹					
MG-C34	4000 PSI	2" OD	1/4" NPT	Chrome ¹					
MG-C63	1/8 LPM*	2" OD	1/4" NPT	Chrome ¹					

MG-C23

OXYGEN REGULATOR REPLACEMENT OUTLET ORIFICE BUSHINGS

For use on Western Flow Gauge Regulators Only

STRAIGHT DIS	SS		STRAIGHT B	ARB
HPG-4949	1/2 LPM	1/4 NPT Male ¹	10-NSP-2	5 LPM 26mm 3/16" Barb to 1/8 NPT Male ¹
10-N80 10-N77	8 LPM #80 15 LPM #77	1/4 NPT Male ¹ 1/4 NPT Male ¹		Daib to 17 o 14F1 Male

^{*}Note: Above orifice configurations match the designated gauge and are required to insure proper performance. When ordering oxygen LPM flow gauges to be used on corresponding Western regulators, the following orifice outlet bushings are required. Choose the proper configuration per chart above.

15 LPM 8 LPM 1/2 LPM 10-N77 10-N80 HPG-4949 10-N77-2

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CYLINDER CONNECTIONS 3000 PSI max				
NUTS	CGA-540			
M54-1	Oxygen Cylinder Nut ¹			
M54-10	Handtight Nut (Brass) ¹	∫ M54-1	M54-10P	
M54-10P	Handtight Nut w/ Plastic Grip ¹	A P		

NIPPLES CGA-540			
M54-2	1/4" NPT Male x 2" ¹		
M54-3	1/4" NPT Male x 2" w/Filter ¹		
M54-4	1/4" NPT Male x 2-1/2" ¹		
M54-5	1/4" NPT Male x 2-1/2" w/Filter ¹		
M54-5CV	1/4" NPT Male x 2-1/2" Check and Filter ¹		M54-5
M54-6	1/4" NPT Male x 3" ¹		
M54-7	1/4" NPT Male x 3" w/Filter ¹		
M54-11	1/4" NPT Male x 2-1/2" w/Soft Tip ¹		M54-11
M54-12	Soft Tip Replacement		M34-11
M54-13	1/4" NPT Male x 4-1/2" ¹		

COUPLER TEES

w/Cylinder Valve Outlets 3000 PSI Max

M54-16 CGA-540 w/Female Inlet to 2 Male Outlets¹

M54-17 CGA-540 w/Female Inlet Check Valve to 2 Male Outlets¹



PLUG & CHAIN ASSEMBLY

Provides gas tight seal 3000 PSI Max.

M54-22 CGA-540¹



¹ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

OXYGEN PIN INDEX YOKES & POST VALVES/SPIGOT VALVES 3000 PSI MAX



OXYGEN YOKES	CGA-870 See Page 41 For Yoke Replacement Parts
PART #	OUTLET
M87-1	CGA-540 Male Thread ¹
M87-1H	CGA-540 Handtight Male Thread ¹
M87-4	1/4 NPT Male ¹
M87-4CV	1/4 NPT Male w/Check Valve ¹
M87-4S	1/4 NPT Male w/Swivel ¹
M87-4H	1/4 NPT Male w/Handtight Yoke ¹
M87-8	1/2 NPT Male ¹
M87-8H	1/2 NPT Male w/Handtight Yoke ¹

OXYGEN POST & SPIGOT VALVES

3/4"-16 UNF-2A 1-1/8"-12 UNF-2A



Note: The cylinder valves listed are configured with proper relief mechanisms for use at 2015 and 2216 PSI. Do not use where cylinder pressure ratings exceed these pressures.

POST VALVE SEALS

MPV-4027 Post Valve Seal (Plastic Clip Dust Cover) See Photo Above¹

TEST GAUGE OXYGEN



¹ ▲ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

DISS NUTS, NIPPLES AND BODY ADAPTORS

	NIPPLES	BODY ADAPTORS
AIR DISS 1160A 3/4-16 Th	read 200 PSI Maximum	
7/8" Hex¹ Handtight (Brass)¹ Handtight/Yellow Collar¹ Handtight/Black/White Collar¹ M16-2P M16-1 M16-2PC	M16-8 M16-9 1/8" NPT Male w/Wrench Flats and O-Ri M16-10 1/4" NPT Male w/Wrench Flats and O-Ri M16-11 1/4" NPT Male w/Wrench Flats and O-Ri M16-12 1/4" Hose Barb 1 1/4" Hose Barb w/O-Ring 1 M16-14 3/16" Hose Barb w/O-Ring 1 M16-15 3/16" Hose Barb w/O-Ring 1 O-Ring Replacement M16-9 M16-12	M16-6 1/8" NPT Female ¹
L DICC 11004 0/41/7	L cooperation :	
7/8" Hex ¹ Handtight (Brass) ¹ Handtight/Black Collar ¹ M12-2 M12-2P	M12-8 1/8" NPT Male w/Wrench Flats 1/8" NPT Male w/Wrench Flats and O-Ri M12-10 1/4" NPT Male w/Wrench Flats 1/4" NPT Male w/Wrench Flats and O-Ri M12-11 1/4" NPT Male w/Wrench Flats and O-Ri M12-12 1/4" Hose Barb 1/4" Hose Barb M12-13 3/16" Hose Barb M12-14 O-Ring Replacement M12-11	
OXIDE DISS 1040A 3/4-	16 Thread 200 PSI Maximum	
7/8" Hex ¹ Handtight (Brass) ¹ Handtight/Blue Collar ¹ M04-2 M04-2P	M04-8 1/8" NPT Male w/Wrench Flats 1/8" NPT Male w/Wrench Flats and O-R M04-10 1/4" NPT Male w/Wrench Flats 1/4" NPT Male w/Wrench Flats and O-R M04-11 1/4" NPT Male w/Wrench Flats and O-R M04-12 1/4" Hose Barb 1/4" Hose Barb M04-13 1/4" Hose Barb w/O-Ring M04-14 3/16" Hose Barb w/O-Ring M04-15 M04-16 O-Ring Replacement Retaining Ring M04-13	M04-6 1/8" NPT Female ¹
	7/8" Hex¹ Handtight (Brass)¹ Handtight/Yellow Collar¹ Handtight/Black/White Collar¹ M16-2P M16-1 M16-2PC DISS 1120A 3/4-16 The 7/8" Hex¹ Handtight (Brass)¹ Handtight/Black Collar¹ M12-2 M12-2P OXIDE DISS 1040A 3/4-17/8" Hex¹ Handtight (Brass)¹ Handtight (Brass)¹ Handtight (Brass)¹ Handtight (Brass)¹ Handtight (Brass)¹ Handtight (Brass)¹ Handtight/Blue Collar¹	MIC-2PC MIC-12 MIC-12 MIC-12 MIC-14 MIC-12 MIC-14 MIC-12 MIC-14 MIC-15 MIC-16 MIC-16 MIC-16 MIC-16 MIC-16 MIC-17 MIC-17 MIC-17 MIC-18 MIC-18 MIC-18 MIC-19 MIC-18 MIC-19 MIC-18 MIC-19 MIC-18 MIC-19 MIC-19

Important: Hose used on straight barbed hose fittings should be secured with crimp-style hose ferrules or proper clamps. Nipples with O-Ring should be used when using handlight nuts to insure proper seal.

¹ **_____WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

DISS NUTS, NIPPLES AND BODY ADAPTORS

NUTS	715, NIPPLES AN	NIPPLES		ΔΓ	DAPTORS
SUCTION	DISS 1220 3/4-16 T				DAI 10K5
M22-1	7/8" Hex ¹	M22-8	1/8" NPT Male Broached for A	Allen Wrench 1	22-4 1/8" NPT Male ¹
M22-1 M22-2	Handtight (Brass) ¹	M22-8 M22-9	1/8" NPT Male Broached w/C		22-5 1/4" NPT Male ¹
M22-2 M22-2P	Handtight/White Coll		1/4" NPT Male w/Wrench Fla	. •	22-6 1/8" NPT Female ¹
M22-2P M22-2PC	Handtight/Yellow Coll		1/4" NPT Male w/Wrench Fla		22-7 1/4" Hose Barb ¹
M22-2PC	nanatignt/ tellow Coll	M22-11	1/4" Hose Barb ¹	is a O-king.	22-7 1/4 Hose barb
	1 3	M22-12 M22-13	1/4" Hose Barb w/O-Ring ¹		M22-5
	M22-2P	M22-16	O-Ring Replacement		MZZ 3
		M22-18	5/16" Hose Barb ¹		
		M22-19	5/16" Hose Barb O-Ring ¹		
İ	M22-2	M22-19	_		
	-	الرباك	M22-9	M22-12	
EVACUATIO	N DISS 2220 7/8-	14 Thread			
M26-1	1" Hex ¹	M26-18	5/16" Hose Barb ¹	M	26-5 1/4" NPT Male ¹
	1 110%	M26-19	5/16" Hose Barb w/O-Ring ¹		
	0	M00-5	Retaining Ring ²		
	M26-1		- Part		M26-5
			M26-19	M00-5	
OXYGEN/C	CARBON DIOXIDE MIX	TURES Carbon D	Dioxide <7% DISS 1200A	3/4-16 Thread	l 200 PSI Maximum
M20-1	7/8" Hex ¹	M20-12	1/4" Hose Barb ¹		20-5 1/4" NPT Male ¹
		M20-13	1/4" Hose Barb w/O-Ring ¹		
	M20-1	M20-16	O-Ring Replacement	M20-13	M20-5
		DE-OXYGEN MIX			4-16 Thread 200 PSI Maximum
M08-1	7/8" Hex ¹	M08-8	1/8" NPT Male w/Wrench Fla		08-5 1/4" NPT Male ¹
M08-2	Handtight (Brass) ¹	M08-9	1/8" NPT Male w/Wrench Fla	ts and O-Ring ¹	
	-	M08-12	1/4" Hose Barb ¹		M08-5
	M08-1	M08-13	1/4" Hose Barb w/O-Ring ¹		THE STATE OF THE S
		M08-16	O-Ring Replacement	M08-12	
HELIUM HI	LIUM OXYGEN MIXTU	RFS Helium > 20	0% DISS 1060A 3/4-16 T	hread 200 PS	il Maximum
M06-1	7/8" Hex ¹	M08-16	O-Ring Replacement		06-5 1/4" NPT Male ¹
M06-2	Handtight (Brass) 1				
	M06-1				M06-5
	YGEN MIXTURES - HE				1 (41) 107 1 1 1
M18-1	7/8" HEX ¹	M18-12	1/4" Hose Barb ¹	M	18-5 1/4" NPT Male ¹
M18-2	Handtight (Brass) 1	M18-13	1/4" Hose Barb w/O-Ring ¹		
SPECIAL MI		3/4-16 Thread	200 PSI Maximum		00 F 1 / 411 NIDT 14 1 1
M02-1	7/8" Hex ¹	M02-12	1/4" Hose Barb ¹	M	02-5 1/4" NPT Male ¹
M02-2	Handtight (Brass) ¹	M02-13	1/4" Hose Barb w/O-Ring ¹		
	H ====	M02-16	O-Ring Replacement		
	M02-1	M00-5	Retaining Ring ²	M00-5	M02-5
l			M02-12	71100 3	_
			Manager and the second		

Important: Hose used on straight barbed hose fittings should be secured with crimp-style hose ferrules or proper clamps. Nipples with O-Ring should be used when using handtight nuts to insure proper seal.

38

¹ **______ WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

NIPPLES		EXTENS	IONS .	Y-CON	NECTOR
3000 PSI A	Max	3000 PS	Max	1000 P	SI Max
MA-26	1/8" NPT Male to $1/8$ " NPT Male x $1-1/2$ "	MA-32	1/8" NPT Female to 1/8" NPT Male x 2-1/2"	MY-5	
MA-27	1/8" NPT Male to $1/8$ " NPT Male x $2-1/2$ "	MA-33	1/4" NPT Female to $1/4$ " NPT Male x 3 "		1/4" NPT Female Outlets ¹
MA-28	1/8" NPT Male to 1/8" NPT Male x $3-1/2$ "			MY-6	1/8" NPT Female Inlet to
MA-29	1/4" NPT Male to 1/4" NPT Male x 1-1/2" 1				1/4" NPT Female Outlets 1
MA-30	1/4" NPT Male to 1/4" NPT Male x 2-1/2" 1		A-		
4	MA-30		MA-33		MY-6

CROSSES	S & TEES 1000 PSI Max Chrome-Plat	ed Brass		
BRANCH	I TEES	STREET	TEES	CROSSES
MA-18	1/8" NPT (3) Female ¹	MA-22	1/8" NPT Female (2) x 1/8" NPT Male ¹	MA-24 1/8" NPT Female x 4 ¹
MA-19	1/4" NPT (3) Female ¹	MA-23	1/4" NPT Female (2) x 1/4" NPT Male ¹	MA-25 1/4" NPT Female x 4 ¹
MA-20	1/8" NPT Male to 1/8" NPT (2) Female ¹			
MA-21	1/4" NPT Male to 1/4" NPT (2) Female ¹			
	MA-19		MA-23	MA-25

ADAPTOR BLOCKS Chrome-Plated Brass 90° w/PIPE THREADS TEE w/PIPE THREADS POWER TAKE-OFF ADAPTOR							
3000 PSI Max	3000 PSI Max	200 PSI Max					
MGA-2 1/8" NPT Female to 1/4" NPT Female ¹ MGA-4 1/4" NPT Female to 1/4" NPT Female ¹ MGA-5 1/4" NPT Male to 1/4" NPT (2) Female ¹	MT-2 1/4" NPT (3) Female x 2" MT-3 1/4" NPT (3) Female x 3-1/8"	MT-5 1/8" NPT Male to 1/8" NPT (2) Female ¹ M24-60 Demand Valve Assembly (MT-5 & M24-40) ¹					
MGA-5	☐ () (MT-2	MT-5 M24-60					

BUSHIN	GS, COUPLERS, PLUGS 3000 PSI N	ax.	
BUSHIN	GS: PIPE x PIPE THREAD	COUPLERS: PIPE x PIPE THREAD	PLUGS: PIPE THREAD
MA-1	1/8" NPT Male to 1/8" NPT Male ¹	MA-4 1/8" NPT Female to 1/8" NPT Female ¹	MA-35 7/16" Hex to 1/8" NPT Male ¹
MA-2	1/4" NPT Male to 1/4" NPT Male ¹	MA-5 1/4" NPT Female to 1/4" NPT Female ¹	MA-36 5/8" Hex to 1/4" NPT Male ¹
MA-3	1/4" NPT Male to 1/8" NPT Male ¹	MA-6 1/4" NPT Female to 1/8" NPT Female ¹	
MA-7	1/8" NPT Female to 1/8" NPT Male ¹		AMERICANIA
MA-8	1/4" NPT Female to 1/4" NPT Male ¹	MA-5	MA-36
MA-9	1/4" NPT Female to 1/8" NPT Male ¹	MA 5	
MA-10	1/8" NPT Female to 1/4" NPT Male ¹	ANGLE ADAPTORS 1000 PSI Max.	
MA-37	1/8" NPT Male to 1/8" NPT Male ¹	90° PIPE x PIPE THREAD	45° PIPE x PIPE THREAD
	MA-3		MA-16 1/8" NPT Female to 1/8 NPT Male MA-17 1/4" NPT Female to 1/4 NPT Male
ANGLE A	ADAPTORS 1000 PSI Max.		
90° PIP	E × PIPE THREAD		MA-16
MA-11	1/8" NPT Female to 1/8" NPT Female ¹ MA-11	MA-14	WHITE
MA-12	1/4" NPT Female to 1/4" NPT Female ¹		

^{1 ▲} WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

NPT MAL	E TO BARB	NPT FEM	ALE TO BARB	HOSE B	ARB SPLICERS
MH-5 MH-6 MH-7 MH-8 MH-9 MH-10	1/8" NPT Male to 1/8" Hose Barb ¹ 1/8" NPT Male to 3/16" Hose Barb ¹ 1/8" NPT Male to 1/4" Hose Barb ¹ 1/4" NPT Male to 1/8" Hose Barb ¹ 1/4" NPT Male to 3/16" Hose Barb ¹ 1/4" NPT Male to 1/4" Hose Barb ¹	MH-12 MH-13 MH-14 MH-15 MH-16 MH-17	1/8" NPT Female to 1/4" Hose Barb ¹ 1/4" NPT Female to 1/8" Hose Barb ¹ 1/4" NPT Female to 3/16" Hose Barb ¹ 1/4" NPT Female to 1/4" Hose Barb ¹	MH-2 MH-1	3/16" Hose Barb Splicer ¹ 1/4" Hose Barb Splicer ¹ MH-1
MH-11 MH-20	1/4" NPT Male to 3/8" Hose Barb ¹ 1/4" NPT Male to 5/16" Hose Barb ¹ MH-10	MH-18 MH-21	1/4" NPT Female to 3/8" Hose Barb ¹ 1/4" NPT Female to 5/16" Hose Barb ¹ MH-17	F	MH-2
	MH-8				

HOSE FERRULES AND CRIMPERS

FERRULE	FITS WEST	TERN HOSE	CRIMPER
M-624-A	MH4-XXN	1/4" ID Non-Conductive	CR-100 ¹
M-625-A	MH4-XXA	1/4" ID Conductive	CR-100 ¹
M-625-A	MH4-XXC	1/4" ID Non-UL	CR-100 or MC-4 ¹
M-620	MH3-00N	3/16" ID Non-Conductive	MC-4 ¹
M-K-K	MH5-22C	5/16" ID Non-UL	CR-100 or MC-4 ¹
M-K-K	MH5-26N	5/16" ID Non-Conductive	CR-100 or MC-4 ¹

(See pages 48-49 for Custom Hose Assemblies and Hose Identification.)



Nickel Plated Brass Ferrules, for use with 200 PSIg Max. Hose assembly.

FERRULES 1

MODEL	ID	LENGTH	CRIMPER*
M-620	.410	1/2"	C-5A
M-769	.450	9/16"	C-5A and CR-100
M-4750	.478	11/16"	CR-100
M-7322	.500	1"	CR-100
M-624-A	.500	5/8"	CR-100
M-624	.500	1/2"	CR-100
M-7323	.525	1"	CR-100
M-625-A	.525	5/8"	MC-4
M-625	.525	1/2"	MC-4
M-626	.548	31/64"	MC-4
M-K-K	.575	31/64"	MC-4

^{*}Recommended Crimper selection based on Ferrule outside diameter (OD). Hose type and size may affect Ferrule and Crimper selection. All hose assemblies should be tested prior to using.

CRIMPERS

	MODEL	DESCRIPTION
Vise-Grip® Style	MC-4 ²	For 3/16" and 1/4" ID Hose

¹ **____WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

² WARNING: This product can expose you to chromium (hexavalent compounds), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

MEDICAL AIR - NITROGEN - NITROUS OXIDE - SUCTION

NON-COF	RROSIVE GASES 1/4 N	PT THREADS	5 200 PSI Max Inlet		
DEMAND	VALVES	CHECK VAI	VES	REVERS	E FLOW VALVES
MEDICAL	AIR	MEDICAL A	IR	NON-C	ORROSIVE GASES
M16-40	Flow @ 50 PSI 1/8" NPT Male ¹ 1/4" NPT Male ¹	50 PSI M16-4CV M16-5CV	DISS 1160-A X 1/8" NPT Male ¹ DISS 1160-A 1/4" NPT Male ¹ M16-4CV	MCV-4	Flow into 1/4" NPT male Inlet to 1/4" NPT male MCV-4
	N I Flow @ 175 PSI 1/4" NPT Male ¹	NITROGEN 50 PSI M12-5CV	-47	MCV-5	Flow into 1/4" NPT Male Inlet to 1/4" NPT Female Outlet ¹
M04-40	M12-50 OXIDE Flow @ 50 PSI 1/8" NPT Male ¹ 1/4" NPT Male ¹			MCV-6	Flow into 1/4" NPT Female to 1/4" NPT Male Outlet ¹
SUCTION 180 SCFH	Flow @ 50 PSI 1/4" NPT Male ¹			MCV-7	Flow into 1/4" NPT Female inlet to 1/4" NPT Female Outlet 1
WAGD M26-50	1/4" NPT Male ¹			_	MCV-7
CARBON					
M08-50	1/4" NPT Male ¹				
MD-1 MD-2 MD-3	1/8" NPT Male to 1/8" NF 1/4" NPT Male to 1/8" NF 1/4" NPT Male to 1/4" NF	PT Female ¹		o clear forwa	alve when plunger hole depth ard travel of plunger when check for full flow after installation.

DEMAND VALVES

Gas is checked w/out connection being made. Flow occurs when proper connection is made.

CHECK VALVES

Gas flows into and through the DISS portion of the fitting and is checked in the opposite direction. Does not operate as a demand valve.

REVERSE CHECK VALVES

For inline gas flow. Prevents backflow of gases.

¹ **WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

PIN INDEX YOKES 3000 PSI MAX



YOKES

IOILES					
Yokes rated 3000	PSI Max				
Nitrogen	M960-4	1/4 NPT Male ¹	Nitrous Oxide	M91-1	CGA-326 Male ¹
CGA-960	M960-8	1/2 NPT Male ¹	CGA-910	M91-4	1/4 NPT Male ¹
				M91-4CV	1/4 NPT Male w/Check Valve ¹
				M91-8	1/2 NPT Male ¹
Nitrous Oxide &	Oxygen mixture	s	Carbon Dioxide & (Oxygen mixture	s
N ₂ O 47.5% to 52.5	%		CO ₂ Not over 7%	M88-4	1/4 NPT Male ¹
CGA-965	M965-4	1/4 NPT Male ¹	CGA-880		
	M965-8	1/2 NPT Male ¹			
Carbon Dioxide 8	Carbon Dioxide	& Oxygen mixtures	Helium & Oxygen ı	mixtures	
CO ₂ over 7%	M94-4	1/4 NPT Male ¹	He Not over 80%	M89-4	1/4 NPT Male ¹
CGA-940	M94-8	1/2 NPT Male ¹	CGA-890		
Helium & Helium	& Oxygen mixto	ures	Non-corrosive, Non	n-flammable mix	ctures
He Not over 80%	M93-4	1/4 NPT Male ¹	CGA-973	M973-1	CGA-500 Female ¹
CGA-930				M973-4	1/4 NPT Male ¹
Medical Air					
CGA-950	M95-1	CGA-346 Male ¹			
	M95-4	1/4 NPT Male ¹			

See CGA chart on inside back cover of catalog.

YOKE REPLACEMENT PARTS







M86-4CV-1 1/4 NPT Male Stem w/Check Valve¹



M86-4-1A 1/4 NPT Male Stem¹





M87-1-1A CGA-540 Stem¹



MCP370 Handtight Yoke Screw¹



Swivel Yoke Body Only – Oxygen



2544-2 Brass and Viton Washer¹



MW-87 Halar Washer

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

M34-22


MEDICAL AIR NITROUS OXIDE NITROGEN HELIUM

MEDICAL A		CO3 CXIDI	E MITROGEM TIELIOM		
COUPLER TEES	Nitrous Oxide	Medical Air	Nitrogen or Helium		
w/Cylinder Val	ve Inlet and Two	Same CGA Outl	ets 3000 PSI Max		
M32-16	CGA-326	Nitrous Oxide w	v/FemaleInlet to 2 Male Outlets ¹		
M34-16	CGA-346	Medical Air w/I	Female Inlet to 2 Male Outlets ¹		
M58-16	CGA-580	Nitrogen or Heli	ium w/Male Inlet to 2 Female Outlets ¹		
PLUG & CHAIN	ASSEMBLY Med	dical Air			M34-16
Provides gas tight	seal	3000 PSI Max.			
M34-22	CGA-346	Medical Air ¹		~	

^{1 ▲} WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CYLINDER CONI	NECTIONS 3000 PSI Max		
NUTS		NIPPLES	
MEDICAL AIR	CGA-346	M34-4	1/4" NPT Male, 2-1/2" Length ¹
830-14 NGO RH	Internal Thread	M34-5	1/4" NPT Male, 2-1/2" Length w/Filter ¹
M34-1	1-1/8" Hex Nut ¹	M34-5CV	1/4" NPT Male, 2-1/2" Length w/Filter and Check Valve ¹
M34-10	Handtight (Brass) ¹	M34-11	1/4" NPT Male, 2-1/2" Length w/Soft Tip ¹
M34-10P	Handtight w/Black Plastic Grip ¹	M34-12	Soft Tip Replacement
	M34-10	M34-10P	M34-5
NITROGEN	OR HELIUM CGA-580	M58-3	1/4" NPT Male, 2-1/2" Length w/Filter ¹
960-14 NGO RH	External Thread	M58-3CV	1/4" NPT Male, 2-1/2" Length w/Filter and Check Valve ¹
M58-1	1-1/8" Hex Nut ¹	M58-5	1/4" NPT Male, 3" Length w/Filter ¹
M58-10	Handtight (Brass) 1	M58-7	1/4" NPT Male, 3-1/2" Length w/Filter ¹
M58-10P	Handtight w/Black Plastic Grip ¹	M58-11	1/4" NPT Male, 3" Length w/Filter w/O-Ring ¹
		M58-13	1/4" NPT Male, 3" Length w/Soft Tip ¹
	M58-1		M58-5
NITROUS OXIDI	CGA-326	M32-4	1/4" NPT Male, 2-1/2" Length ¹
830-14 NGO RH		M32-5	1/4" NPT Male, 2-1/2" Length w/Filter ¹
M32-1	1-1/8" Hex Nut ¹	M32-5CV	1/4" NPT Male, 2-1/2" Length w/Filter and Check Valve ¹
M32-10	Handtight (Brass) ¹	M32-11	1/4" NPT Male, 2-1/2" Length w/Soft Tip ¹
M32-10P	Handtight w/Black Plastic Grip ¹		,
	M32-1		M32-11
CARBON DIOXI	DE CGA-320	M30-2	1/4" NPT Male, 2" length w/Washer ¹
830-14 NGO RH	Internal Thread	M30-4	1/4" NPT Male, 2-1/2" length w/Washer ¹
M30-1	1-1/8" Hex Nut ¹	M30-5	1/4" NPT Male, 2-1/2" length w/Washer and Filter ¹
M30-10P	Handtight w/Black Plastic Grip ¹	M30-5CV	1/4" NPT Male, 2-1/2" length w/Washer/Filter/Check Valve ¹
	M30-1		M30-4
MEDICAL BREAT	HING MIXTURES CGA-280	M28-4	1/4" NPT Male, 2-1/2" Length ¹
750-14 NGO RH	Internal Thread	M28-5	1/4" NPT Male, 2-1/2" Length w/Filter ¹
M28-1	1-1/8" Hex Nut ¹	M28-11	1/4" NPT Male, 2-1/2" Length w/Soft Tip ¹
M28-10P	Handtight w/Black Plastic Grip ¹	M28-12	Soft Tip Replacement
	M28-1		M28-5
MEDICAL MIXTU	JRES CGA-500	M50-3	1/4" NPT Male, 2-1/2" Length w/Filter ¹
.830-14 NGO RH	External Thread	M50-7	1/4" NPT Male, 3-1/2" Length w/Filter ¹
M50-1	1-1/8" Hex Nut ¹	M58-11	1/4" NPT Male, 3" Length w/O-Ring ¹
M50-10	Handtight (Brass) 1		•
M50-10P	Handtight w/Black Plastic Grip ¹		
	M50-1		M50-3

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

INLET CONNECTIONS	OXYGEN DISS 1240	MEDICAL AIR DISS 1160-A	SUCTION DISS 1220
1/8" NPT Female Without Connector ¹	M24-70	M16-70	M22-70
1/4" NPT Female Without Connector ¹	M24-71	M16-71	M22-71
DISS Hex Nut and Nipple ¹	M24-74	M16-74	M22-74
DISS Handtight Nut and Nipple*1	M24-75	M16-75	M22-75
Chemetron® Adaptor*1	M24-76	M16-76	M22-76
Ohmeda® Adaptor* ¹	M24-77	M16-77	M22-77
Oxequip® Adaptor ¹	M24-78	M16-78	M22-78
Puritan® Adaptor ¹	M24-79	M16-79	M22-79







Example:

M24-75 2 outlet style with DISS Handtight Nut & Nipple Connection

⁵⁰ PSI Nominal Working Pressure

^{*}Canadian color-coding available on most models, add "C" to part number. Example: M24-77C Product not to scale.

¹▲WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Ch	emetron [®]	OR Y-Coupl	Ohmeda® OR Y-Couplers					
INLET CONNECTIONS	INLET CONNECTIONS OXYGEN			SUCTION	OXYGEN	MEDICAL AIR	NITROUS OXIDE	SUCTION
1/8" NPT Female Without Connector ¹	OR101-70	OR102-70	OR103-70	OR104-70	OR201-70	OR202-70	OR203-70	OR204-70
DISS Hex Nut and Nipple ¹	OR101-74	OR102-74	OR103-74	OR104-74	OR201-74	OR202-74	OR203-74	OR204-74
DISS Handtight Nut and Nipple*1	OR101-75	OR102-75	OR103-75	OR104-75	OR201-75	OR202-75	OR203-75	OR204-75
Chemetron® Adaptor*	OR101-76	OR102-76	OR103-76	OR104-76	OR201-76	OR202-76	OR203-76	OR204-76
Ohmeda® Adaptor* ¹	OR101-77	OR102-77	OR103-77	OR104-77	OR201-77	OR202-77	OR203-77	OR204-77

⁵⁰ PSI Nominal Working Pressure *Canadian color-coding available on most models, add "C" to part number. Example: OR201-75C





Example:

OR101-76 Chemetron® OR Y-Coupler with Chemetron® Adaptor



OR201-70



Example:

OR204-76 Ohmeda® OR Y-Coupler with Chemetron® Adaptor

¹ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



1 HOSE TYPE	2 HOSE SIZE	3 GAS SERVICE	4 HOSE COLOR	5 HOSE LENGTH	6 and 7 INLET-OUTLET CONNECTIONS
		1	1		

C = Conductive	4 = 1/4"	Al			Specify Hose Length
N = Non-Conductive	5 = 5/16"	OXY	Oxygen	Designation Default	in Inches
		SUC	Suction	Refer to chart for Canadian Color Codes	Minimum of 12"
		NI	Nitrogen	Canadian Color Codes	
		N ₂	Nitrous Oxide		
		CO	Carbon Dioxide		

HOW TO ORDER:

Western's custom hose assemblies conform to CGA standards. All assemblies are cleaned and poly-bagged for delivery. To order, select your requirements from the above charts.

GAS DESIGNATION/HOSE COLOR CHART

			CANADA				
Gas	Designation		Color	or Designation		Color	
Medical Air	Al	_	Yellow	Al	K	Black	
Oxygen	OX	_	Green	OX	W	White	
Suction	SU	_	White	SU	Υ	Yellow	
Nitrogen	NI	_	Black	NI	_	Black	
Nitrous Oxide	N2	_	Blue	N2	_	Blue	
Carbon Dioxide	CO	_	Grey	CO	_	Grey	

Female Nut and Nipple FNN MBA Male Body Adaptor MP2 Male Pipe 1/8" Female Pipe 1/8" FP2 Ohmeda® Male OHM Puritan® Male PBM Chemetron® Male CHM Blank End BLK Female Handtight **FHT** DBA Demand Body Adaptor Male Pipe 1/4" MP4 FP4 Female Pipe 1/4" OHF Ohmeda® Female MXO Oxequip® Male SHF Schrader® Female Chemetron® Female CHF

Note: For 90° Elbow add "L" before Inlet or Outlet code (LFNN)

HOSE COLOR:

All medical hose color coding will default to USA standards unless otherwise specified.

HOSE LENGTH REQUIREMENTS:

Hose length 12 inches minimum.

Lengths 12" to 60" (1 ft - 5 ft) are available in 6 inch increments only.

Lengths 60" to 300" (5 ft - 25 ft) are available in 12 inch (1 foot) increments only.

Lengths over 300" (25 feet) shall be in 60 inch (5 foot) increments only.

Conductive hose maximum length is 240 inches (20 feet) due to conductivity.

Non-conductive hose maximum length is 600 inches (50 feet).

EXAMPLE

USA Part #: C4OX-120FNNFNN

Canadian Part #: C4OXW120FNNFNN

C Conductive Hose

4 1/4"

OX Oxygen (USA)


OXW Oxygen (Canada)

120 120" Long

FNN Female Nut and Nipple Inlet

FNN Female Nut and Nipple Outlet

INLET-OUTLET CONNECTION OPTIONS





REINFORCED MEDICAL HOSE MH4-24C MH4-16N **BULK HOSE** CONDUCTIVE, UL LISTED 250' REEL **NON-CONDUCTIVE 300' REEL** MH4-04A 1/4" ID, Nitrous Oxide (Blue) MH3-00N 3/16" ID, General (Clear) CONDUCTIVE, UL LISTED MH4-12A 1/4" ID, Nitrogen (Black) MH4-00N 1/4" ID, General (Clear) 250 PSI Working Pressure 1000 PSI Burst Pressure MH4-16A 1/4" ID, Medical Air (Yellow), MH4-04N 1/4" ID, Nitrous Oxide (Blue) Suction for Canada MH4-12N 1/4" ID, Nitrogen (Black) CONDUCTIVE, NON-UL LISTED MH4-22A 1/4" ID, Suction (White), MH4-16N 1/4" ID, Medical Air (Yellow), Oxygen for Canada 400 PSI Working Pressure Suction for Canada 1600 PSI Burst Pressure MH4-24A 1/4" ID, Oxygen (Green) MH4-22N 1/4" ID, Suction (White), Oxygen for Canada **CONDUCTIVE, NON-UL LISTED 300' REEL** NON-CONDUCTIVE MH4-24N 1/4" ID, Oxygen (Green) 250 PSI Working Pressure MH4-04C 1/4" ID, Nitrous Oxide (Blue) 1000 PSI Burst Pressure MH4-08C 1/4" ID, Carbon Dioxide (Grey) MH4-12C 1/4" ID, Nitrogen (Black) MH4-16C 1/4" ID, Medical Air (Yellow), Suction for Canada 1/4" ID, Suction (White), Oxygen MH4-22C for Canada MH4-24C 1/4" ID, Oxygen (Green) MH5-22C 5/16" ID, Suction (White), Oxygen for Canada

¹ ▲ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Standard. Medical Cylinder Pin Index Valve Yokes And Threaded Cylinder Connections for Pressures up to 3000 PSI (20680 kPa)

DISS Connections 200 PSI Max

Cylllide	Conneci	ons for Pressures up to 3000 PSI (20080 KPa)	DISS Connections 200 PSI Max.				
YOKE	THREAD	GAS ASSIGNMENT	DISS	GAS ASSIGNMENT			
870	540	Oxygen	1020A	Special Mixtures (for limited experimental applications)			
880	280	Carbon Dioxide & Oxygen Mixture (CO ₂ not over 7.0%)	1040A	Nitrous Oxide			
890	280	Helium & Oxygen Mixtures (Helium not over 80%)	1060A	Helium, Helium-Oxygen Mixtures (Helium over 80%)			
890	280	Nitrogen & Oxygen Mixtures (O ₂ over 23.5%)	1080A	Carbon Dioxide, Carbon Dioxide-Oxygen Mixtures (Carbon Dioxide over 7%)			
890	280	Xenon & Oxygen Mixtures (O ₂ over 20%)	1100A				
910	326	Nitrous Oxide	1120A	Nitrogen			
930	500	Helium & Helium Oxygen Mixtures (Helium over 80%)	1140A	Ethylene			
940	500	Carbon Dioxide & Carbon Dioxide Oxygen Mixtures (CO ₂ over 7%)	1160A	Medical Air			
950	346	Medical Air	1180A	Oxygen-Helium Mixtures			
				(Helium not over 80%)			
965	280	Nitrous Oxide & Oxygen Mixture (N ₂ O 47.5 to 52.5%)	1200A	Oxygen-Carbon Dioxide Mixtures (Carbon Dioxide not over 7%)			
960	580	Nitrogen	1220	Suction (Vacuum)			
973	500	Nonflammable, Noncorrosive Gas Mixtures Labeled as	1240	Oxygen			
		Drugs or Medical Devices and not having another Connection Assignment *	2020	Nitrous Oxide-Oxygen Mixtures (Nitrous Oxide 47.5% to 52.5%)			
		*This is the yoke connection for gas mixtures used with the CGA 500, which are not otherwise assigned.	2040	Oxygen or Oxygen Enriched Air Mixtures			
973	500	Carbon Dioxide, Oxygen, Nitrogen	2060	Xenon			
973	500	Clinical Blood Gas Mixtures	2080	Instrument Air			
973	500	Lung Diffusion Mixtures	2220	Gas Evacuation			

MEDICAL COLOR CODES

Carbon Dioxide Nitrous Oxide/Oxygen Mix Carbon Dioxide/Oxygen Mix Helium Helium/Oxygen Mix Medical Air Nitrogen Nitrogen/Oxygen Mix Nitrous Oxide Xenon/Oxygen Mix Oxygen

EASY GUIDE TO WESTERN MEDICA PART NUMBERS

MEDICAL GAS FITTINGS

HOW IT WORKS:

M24-1 M54-1 M87-1

(M=Medical 24=DISS 1240 1= Part#) (M=Medical 54=CGA-540 1=Part #) (M=Medical 87=CGA-870 1=Part #) Medical Oxygen DISS Nut Medical Oxygen Cylinder Nut Medical Oxygen Yoke

MEDICAL GAS REGULATORS

Part Number- M1-540-FG

(M=Medical 1=Single Stage 540=CGA-540 FG=Flow Gauge) Description: Medical Single-Stage Oxygen Flow Gauge Regulator

WESTERN MAKES IT EASIER...

Interested in the benefits of gas manifolds, but not interested in the complexities?

Let Western's team of experienced sales, technical support, and engineering be your partner.

Western knows the gas manifold specification process is complex. The Western Manifold Configurator helps to demystify the complexity of selecting the right manifold for your application. The Western Manifold Configurator will simplify the complex.



Visit www.westernenterprises.com/configurator/ to activate our team that will work alongside you to determine your manifold needs, recommend a specific solution, and deliver a quality product to achieve success during the process.

LIAMITED	VA/A DD	ARITY	INFORM	ATION
	VVARR		HALOKIA	

Western Enterprises warrants components of the products listed in this catalog to be free of manufacturing defects for the time periods listed below. Please read the Limited Warranty information included with the product for complete details.

Valve Integrated Pressure Regulators	5 years (OxyTote, OxyTote NG, AirTote)
Click Style Regulators	2 years (OPA, OPL, Titan)
MLR Regulators	1 year
Medical Regulators/Flowgauges	6 years (M1, M2)
Flowmeters	6 years (FM, FMA)
Click Style Flowmeters	2 years (CFM)
Transfill Equipment	1 year
Hoses	5 years
Fittings/Quick Connects	90 days
Carts, Racks, Stands	90 days
Bags	90 days

Part /	D (.)	Part /	D (.)	Part /	D (.)	Part /	D (.)	Part /	D (.)	Part /	D (.)	Part /	D (.)	Part /
	Page(s)	Model #	Page(s)	Model #	Page(s)	Model #	Page(s)	Model #	Page(s)	Model #	Page(s)	Model #	Page(s)	Model # Page(s
10-N77		M12-2P		M1-870-PED		M24-30		M54-10		MA-27		1	49	MQ3012
10-N80		M12-5		M1-870-PG		M24-31		M54-10P		MA-28			49	MQ3022
10-NSP-2		M12-50		M1-870-PGE		M24-32		M54-11		MA-29		MH4-04A .		MQ3032
2544-2		M12-5CV		M1-880-PG		M24-33		M54-12		MA-3		MH4-04C.		MQ3043
C4OX-240FHT		M12-7		M1-890-P		M24-34		M54-13		MA-30		MH4-04N.		MQ4012
M101		M12-8		M1-890-PG		M24-35		M54-16		MA-32		MH4-08C.		MQ4022
M102		M1-280-P		M1-910-P		M24-36		M54-17		MA-33		MH4-12A.		MQ4032
M103		M1-280-PG.		M1-910-PG		M24-37		M54-2		MA-35		MH4-12C.		MQ4043
M104		M1-280-PGE		M1-930-PG		M24-39		M54-22		MA-36 MA-37		MH4-12N.		MR-870-15FG1
-M106 -M107		M12-9		M1-940-12F		M24-4		M54-3 M54-4		MA-3/		MH4-16A .		MR-870-P1
		M1-320-12F		M1-940-P		M24-40						1	49	MR-870-PED1
-M108 -M109		M1-320-12F		M1-940-PG		M24-41		M54-5		MA-5			49	MSG1805001
		M1-320-P		M1-950-15F		M24-42		M54-5CV M54-6		MA-6		MH4-22A.		MSG1805801
M201		M1-320-PG.		M1-950-P		M24-43				MA-7		MH4-22C.		MSH152801
M202		M1-320-PGE		M1-950-PG		M24-45P		M54-7		MA-8 MA-9		MH4-22N.		MSH153201
M203		M1-320-PGF		M1-950-PGI M1-960-P		M24-45PGF		M58-1		MC-4		MH4-24A .		MSH153261
M204 M208		M1-326-P		M1-960-P		M24-45PW		M58-10 M58-10P				MH4-24C. MH4-24N.	49	MSH153461 MSH155001
-M206		M1-326-PG. M1-326-PGF		M1-900-FG		M24-45PY . M24-47		M58-10F		MCC1		1	49	MSH155401
M603		l .						M58-11				MH5-22C.		
M604		M1-346-15F		M1-973-PGE			32			MCS		1		MSH155801
		M1-346-P		M1A-346-P		M24-48S		M58-16		MCV-1			40	MSH158701
M607		M1-346-PG.		M1A-346-Ph		M24-4CV		M58-3		MCV-2			40	MSH159401
M609		M1-346-PGH		M1A-540-P		M24-5		M58-3CV		MCV-3		MH-8		MSH159601
M802		M1-346-PH.		M1A-540-Ph		M24-50		M58-5		MCV-4		MH-9		MSH159731
M803		M1-500-P		M1A-870-P		M24-5CV		M58-7		MCV-5		MJ-C1		MSH1802801
MR100		M1-500-PG.		M1A-910-P		M24-6		M-620		MCV-6		1	34	MSH1803461
HCC1-M		M1-500-PGE		M1A-940-P		M24-60	39	M-624		MCV-7		MNR-600.	7	MSH1805001
ICC2-M		M1-500-PGH		M1A-950-P		M24-7		M-624-A		MCW-2B		MNR-800.		MSH1805401
HCC3-M		M1-500-PH.		M20-1		M24-70		M-625		MCW-2BC		MNS-602.		MSH1805801
HCC5-M		M1-540-15F		M20-12		M24-71		M-625-A		MCW-3P		MNS-603.		MSH1808701
HCS2-M		M1-540-15F		M20-13		M24-74		M-626		MCW-3PC		MNS-802.		MSH1809501
HCS3-M		M1-540-15F		M20-16		M24-75		M-7322		MD-1		MNS-803.		MSH1809601
HPG-4949		M1-540-8FG	310	M20-5		M24-76		M-7323		MD-2		MPV-4027		MSH4503461
v100-3	42	M1-540-8FM	١13	M22-1	38	M24-77	46	M-769		MD-3		MPV-872		MSH4505401
M00-5	37	M1-540-8FM	۱413	M22-10	38	M24-78	46	M86-4-1A	42	MDM-8	9	MPV-872T.	18	MSH4505801
M00-6	31	M1-540-P	12	M22-11	38	M24-79	46	M86-4CV-1	42	MGA-2	39	MQ101	27	MSH4509601
M00-9	42	M1-540-PED		M22-12	38	M24-8	31	M86-8-1A	42	MGA-4	39	MQ101-1.	27	MSH452801
M02-1	38	M1-540-PED	310	M22-13	38	M24-9	31	M87-1	36	MGA-5	39	MQ101-2.	27	MSH453201
M02-12	38	M1-540-PG.	11	M22-16	38	M2-540-15I	FG10	M87-1-1A	42	MG-C1	34	MQ101-3.	27	MSH453461
M02-13	38	M1-540-PGE	314	M22-18	38	M26-1	38	M87-4	36	MG-C11	45	MQ101-4.	27	MSH455001
M02-16	38	M1-540-PGH	H11	M22-19	38	M26-18	38	M87-4CV	36	MG-C13	45	MQ101-5.	27	MSH455401
M02-2	38	M1-540-PH.	12	M22-2	38	M26-19	38	M87-4S	36	MG-C14	45	MQ101-6.	27	MSH455801
M02-5	38	M1-580A-PH	112	M22-2P	38	M26-5	38	M87-8	36	MG-C15	45	MQ102	28	MSH459401
M04-1	37	M1-580-P	12	M22-2PC	38	M28-1	44	M88-4	42	MG-C16	45	MQ102-1.	28	MSH459501
M04-10	37	M1-580-PG.	11	M22-4	38	M28-10P	44	M89-4	42	MG-C17	45	MQ102-2.	28	MSH802801
M04-11	37	M1-580-PH.	12	M22-5	38	M28-11	44	M91-1	42	MG-C18	45	MQ102-3.	28	MSH803201
M04-12	37	M16-1	37	M22-50	38	M28-12	44	M91-4	42	MG-C19	45	MQ102-4.	28	MSH803261
M04-13	37	M16-10	37	M22-6	38	M28-4	44	M91-4CV	42	MG-C2	34	MQ102-5.	28	MSH803461
M04-14	37	M16-11	37	M22-7	38	M28-5	44	M91-8	42	MG-C23	34	MQ102-6.	28	MSH805001
M04-15	37	M16-12	37	M22-70	46	M30-1	44	M93-4	42	MG-C29	45	MQ104	30	MSH805401
M04-16	37	M16-13	37	M22-71	46	M30-10P	44	M94-4	42	MG-C3	34	MQ104-1.	30	MSH805801
M04-2		M16-14	37	M22-74	46	M30-2	44	M94-8		MG-C30		MQ104-3.	30	MSH809501
M04-2P	37	M16-15	37	M22-75	46	M30-4	44	M95-1	42	MG-C31	34	MQ201	27	MT-23
۸04-4	37	M16-16	37	M22-76	46	M30-5	44	M95-4	42	MG-C32	34	MQ201-1.	27	MT-33
Λ04-40	41	M16-2	37	M22-77	46	M30-5CV	44	M95-8	42	MG-C34	34	MQ201-2.	27	MT-53
Λ04-5		M16-2P		M22-78		M32-1	44	M960-4		MG-C4		MQ201-3.		MTT-1K2
λ04-50		M16-2PC		M22-79		M32-10		M960-8		MG-C5		MQ201-4.		MTT-2K2
Λ04-6		M16-4		M22-8		M32-10P		M965-4		MG-C58		MQ201-5.		MTT-3K2
104-7		M16-40		M22-9		M32-11		M965-8		MG-C6		MQ201-6.		MW-872
۸04-8		M16-4CV		M24-1		M32-16		M973-1		MG-C63		MQ202		MY-13
104-9		M16-5		M24-10		M32-4		M973-4		MG-C67		MQ202-1.		MY-53
106-1		M16-50		M24-11		M32-5		MA-1		MG-C7		MQ202-2.		MY-63
106-2		M16-5CV		M24-12		M32-5CV		MA-10		MG-C75		MQ202-3.		MY-73
ЛО6-5		M16-6		M24-13		M34-1		MA-11		MG-C8			28	NMR-870-15FM1
108-1		M16-7		M24-14		M34-10		MA-12		MG-C9		MQ202-5.		OPA-510
108-12		M16-70		M24-15		M34-10P		MA-13		MGT-540		MQ202-6.		OPA-520
108-13		M16-71		M24-16		M34-11		MA-14		MH-1		MQ203		OPA-540
108-16		M16-74		M24-18		M34-12		MA-15		MH-10		1	29	OPA-810
108-2		M16-75		M24-2		M34-16		MA-16		MH-11		MQ203-2		OPA-820
108-5		M16-76		M24-20		M34-22		MA-17		MH-12		MQ203-3.		OPA-840
108-8		M16-77		M24-24		M34-4		MA-18		MH-13		MQ203-4.		OPA-850
108-9		M16-78		M24-25		M34-5		MA-19		MH-14		MQ203-5.		OPL-821
лоо- 7 Л12-1		M16-79		M24-25		M34-5CV		MA-2		MH-15		MQ203-5.		OR1012
лт2-тл Л12-10		M16-8		M24-27		M-4750		MA-20		MH-16			30	OR101-12
лт2-то Л12-11		M16-9		M24-27		M50-1		MA-20		MH-17		MQ204 MQ204-1.		OR101-12
		M1-870-15F		M24-28 M24-29				MA-21				MQ204-1. MQ204-2.		OR101-22
112-12 112-12		l .				M50-10				MH-18				
Λ12-13		M1-870-15F		M24-2P		M50-10P		MA-23		MH-19		MQ204-3.		OR101-42
		M1-870-5FG M1-870-P		M24-2PA		M50-3		MA-24		MH-2			30	OR101-52
112-14 112-14			12	M24-2PH	31	M50-7	44	MA-25	39	MH-20	40	MQ204-5.	30	OR101-62
۸12-16					2.1	AA5 4 1	2 5	MAA 24	20	사사나 그 1	40	MODOLL	20	OP10170 4
		M1-870-PED		M24-2PW	31	M54-1	35	MA-26	39	MH-21	40	MQ204-6.	30	OR101-704

Part /	
Model #	
OR101-74. OR101-75.	
OR101-76.	
OR101-77.	
OR102 OR102-1	28
OR102-2	
OR102-3	
OR102-4 OR102-5	28 28
OR102-6	28
OR102-70. OR102-74.	
OR102-74.	
OR102-76.	
OR102-77. OR103	
OR103-1	
	29
OR103-3 OR103-4	
OR103-5	29
OR103-6	29
OR103-70. OR103-74.	
OR103-75.	47
OR103-76. OR103-77.	
OR103-77.	
OR104-1	30
OR104-2 OR104-3	30
OR104-3	
OR104-5	30
OR104-6 OR104-70.	
OR104-70.	
OR104-75.	47
OR104-76. OR104-77.	
OR104-77.	
OR201-1	
OR201-2 OR201-3	
OR201-4	27
OR201-5	
OR201-6 OR201-70.	
OR201-74.	47
OR201-75.	
OR201-76. OR201-77.	4/ 47
OR202	
OR202-1	
OR202-2 OR202-3	
OR202-4	28
OR202-5 OR202-6	
OR202-0 OR202-70.	
OR202-74.	47
OR202-75. OR202-76.	
OR202-76.	47
OR203	29
OR203-1 OR203-2	
OR203-2	
OR203-4	
OR203-5 OR203-6	
OR203-70.	
OR203-74.	47
OR203-75 . OR203-76 .	4/ 47
OR203-77.	47
OR204	30
OR204-1 OR204-2	30
,	

Part / Model # Page(s))
•
OR204-330
OR204-7047
OR204-7447
OR204-7547
OR204-7647
OR204-7747
PX-100021
PX-100221
PX-101020
PX-101120
PX-101320
PX-101420
PX-102022
PX-102523
PX-102723
PX-103423
PX-10447
PX-10457
PX-106022
PX-185021
QUE-900136
WPX-5405-119
WPX-5406-119
WPX-8700-1T18
WPX-8700-1W18
WPX-8701-1T19
WPX-8702-1T19
WPX-8703-1T19
WPX-8703-1W19





Medical Products





Global Gas Management Technology

Enterprises | Innovator | **Medica** | Westwinds

875 Bassett Road, Westlake, OH 44145

Toll Free (800)783.7890 / Phone (440)871.2160 / Fax (440)835.8283

ISO 9001/ISO 13485 Registered Quality System

